



autoflex[®]
the high-throughput
MALDI-TOF




BRUKER
DALTONICS[®]

Enabling Life Science Tools Based on Mass Spectrometry™

“The new challenges of the post-genomic era lie before us: recognizing genomic variability, reading the proteomes of hundreds of cell types and their variability, and exploring the function or malfunctions of millions of proteins.”

Meet the challenges of the post-genomic era



HIGH
A THROUGH

We now can read the human genome. We have learned, in principle, how genes “translate” into proteins by the genetic code. But there are millions of “pre-translational modifications (the mutations), of “translational modifications” (the slicing variants), and of “posttranslational modifications” like glycosylation, phosphorylation, cross-linkings and many others. A major challenge!

Reading the human genome was a tool-driven effort: without the development of fast sequencers we could not have celebrated this success. Reading the mutations, the proteomes, the modifications will be another tool-driven effort - and mass spectrometry will play a dominant role.

Demands on the instrumentation for this task are high: automated, high sample throughput mass spectrometers are required, higher mass ranges, and more accurate mass determinations.

With instruments like the autoflex, designed for tens of thousands of MALDI samples per day, Bruker Daltonics will help you meet the challenges of the post-genomic era.

SAMPLE
INPUT



autoflex – trained to be top

The autoflex MALDI-TOF MS is the first MALDI-TOF system thoroughly designed for High Sample Throughput in Genomics and Proteomics.

What does that mean?

Automation

Ultra-fast, automated analysis of thousands of samples, carried out with the high data quality known from the Daltonics FLEX series of MALDI-TOF systems.

Flexibility

Modularity: another characteristic of a Daltonics system. The linear autoflex, tailored for SNP genotyping, can be upgraded with reflector and AutoPSD accessory to make it a powerful protein analysis tool.

Reliability


All system components have been carefully designed and tested for maximum life and minimum wear.

PERFORMANCE




autoflex: a quantum leap
in MALDI-TOF automation
and sample throughput

The autoflex won
the iF Product Design Award 2001
of the German Industrie Forum
Design for its "uncompromising
symbiosis of high-tech and ergonomics".




Flexibility by Design

- Small footprint – saves expensive lab space
- Horizontal, height adjustable probe inlet – adaptable to any process line
- Linear TOF configuration can be upgraded with reflector and PSD equipment




The reflectron increases mass spectrometric resolution



The autoflex adapts to any operating height your lab facilities require

Horizontal feed for sample targets – easy access for either manual or robot driven sample supply



Flexibility in Use: flexControl™

flexControl is Bruker Daltonics' new, universal MALDI-TOF software platform for quick and easy instrument control. The GUI is highly intuitive and provides quick drag-and-drop of information, as well as context-sensitive menus.

flexControl features Bruker Daltonics' proven AutoXecute module for intelligent on-the-fly optimization of acquisition parameters.




AutoXecute methods are pre-selectable and can even cover on-line post-processing, e.g. protein database search.

flexControl gives straight access to all relevant instrument parameters

Seamless Integration

flexControl networks with the Daltonics Proteomics (mapControl™/biotoools™/MASCOT™) and Genomics (genotoools™) software suites.

The "Calibration" tab allows for quick calibration setup



The "General" tab allows to check the system status, including all basic operating parameters

- 1 Online Function Control (OFC)
- 2 Video optics: spots can be manually selected via mouse click!
- 3 Acquisition Control
- 4 Acquisition Parameter Control

As a real high-throughput system, the autoflex features maximum reliability of all system components for many millions of analyses.



Automated Sample Handling

- Safe, automated microtiter plate handling by the Twister

Optimized MALDI Process

- Low-friction xy-table for many tens of millions of speedy “stop and go” operations
- A laser which can provide over 100 million shots without service
- Highly reproducible positioning of the sample spots using our patented AnchorChip technology, which overcomes the need for “spot-hunting”.

The sample tray:
easy-to-handle sample supply



Remote service via ISDN point to point links

Extended Self Diagnostic and Service Capabilities

- Automated self diagnosis routines
- Remote service capability via safe ISDN point-to-point link

Maximum up-time

The built-in remote service capability allows the Bruker Daltonics help desk to access the customer's system and do on-line troubleshooting. Often problems can be solved in this way. In case that hardware faults are localized on-line, the service engineer can take along the spare needed part with him.


reliability

automation

Full Microtiterplate Compatibility

The microtiterplate format of the scoutMTP MALDI targets guarantees for optimum workflow in the lab without compatibility problems.


The Twister Microplate Handler automatically processes up to 80 scoutMTP targets, each carrying up to 1536 samples. Sample targets can be identified via barcode and transponder.



The Twister smoothly interfaces sample preparation in the biochemical lab with the mass spectrometer.

The key to high productivity

The patented AnchorChip coating of the scout MTP targets provides highest sensitivity from smallest sample amounts. Highly reproducible positioning of the sample spots using Anchor-Chips shortens analysis times considerably.



The microtiter plate format of the scout MTP™ MALDI targets guarantees for optimum workflow in the biochemical lab without compatibility problems

1995 – 26 samples

1998 – 384 samples

1999 – 1,536 samples

2000 – 30,000 samples

High sample throughput

Through the last decade, we have trained our FLEX MALDI TOF systems to process a growing number of samples automatically.




Result list in ASCII format

Sample ID	Genotype	Reliability
S_1001	A/T	High
S_1002	A/T	High
S_1003	A/T	High
S_1004	A/T	High
S_1005	A/T	High
S_1006	A/T	High
S_1007	A/T	High
S_1008	A/T	High
S_1009	A/T	High
S_1010	A/T	High
S_1011	A/T	High
S_1012	A/T	High
S_1013	A/T	High
S_1014	A/T	High
S_1015	A/T	High
S_1016	A/T	High
S_1017	A/T	High
S_1018	A/T	High
S_1019	A/T	High
S_1020	A/T	High
S_1021	A/T	High
S_1022	A/T	High
S_1023	A/T	High
S_1024	A/T	High
S_1025	A/T	High
S_1026	A/T	High
S_1027	A/T	High
S_1028	A/T	High
S_1029	A/T	High
S_1030	A/T	High
S_1031	A/T	High
S_1032	A/T	High
S_1033	A/T	High
S_1034	A/T	High
S_1035	A/T	High
S_1036	A/T	High
S_1037	A/T	High
S_1038	A/T	High
S_1039	A/T	High
S_1040	A/T	High
S_1041	A/T	High
S_1042	A/T	High
S_1043	A/T	High
S_1044	A/T	High
S_1045	A/T	High
S_1046	A/T	High
S_1047	A/T	High
S_1048	A/T	High
S_1049	A/T	High
S_1050	A/T	High
S_1051	A/T	High
S_1052	A/T	High
S_1053	A/T	High
S_1054	A/T	High
S_1055	A/T	High
S_1056	A/T	High
S_1057	A/T	High
S_1058	A/T	High
S_1059	A/T	High
S_1060	A/T	High
S_1061	A/T	High
S_1062	A/T	High
S_1063	A/T	High
S_1064	A/T	High
S_1065	A/T	High
S_1066	A/T	High
S_1067	A/T	High
S_1068	A/T	High
S_1069	A/T	High
S_1070	A/T	High
S_1071	A/T	High
S_1072	A/T	High
S_1073	A/T	High
S_1074	A/T	High
S_1075	A/T	High
S_1076	A/T	High
S_1077	A/T	High
S_1078	A/T	High
S_1079	A/T	High
S_1080	A/T	High
S_1081	A/T	High
S_1082	A/T	High
S_1083	A/T	High
S_1084	A/T	High
S_1085	A/T	High
S_1086	A/T	High
S_1087	A/T	High
S_1088	A/T	High
S_1089	A/T	High
S_1090	A/T	High
S_1091	A/T	High
S_1092	A/T	High
S_1093	A/T	High
S_1094	A/T	High
S_1095	A/T	High
S_1096	A/T	High
S_1097	A/T	High
S_1098	A/T	High
S_1099	A/T	High
S_1100	A/T	High

Automated SNP Genotyping

With its high reliability and sample throughput, MALDI-TOF MS is the method of choice for SNP genotyping. The autoflex allows for automated genotyping on an industrial scale within a few seconds per sample.


H63D: C/C homozygous
C282Y: C/C homozygous



Duplex genotyping assay for H63D and C282Y in the HFE gene


The necessary data quality is already achieved in the linear TOF configuration. An intelligent acquisition routine judges the data quality on-the-fly and, if necessary, optimizes the acquisition parameters. The genotools™ software assigns the genotypes and displays the results in an MTP viewer format.

H63D: C/G heterozygous
C282Y: C/T heterozygous



Result overview for 379 samples


Quality control of the results




Automated Protein Identification

Mass spectrometry is the key technology for unlocking the proteome. And in current mass-spectrometric based proteomics strategies, MALDI-TOF MS is the initial screening step – thus the speed and reliability of your MALDI-TOF system will decide on your high throughput success in protein identification!

Protein identification as batch process: the biotools software automatically classifies the quality of the search result for each sample spot and displays a color coded result overview.




High quality (most important: sensitivity) of MALDI-TOF data is crucial for high success rates in Protein ID. Here: 300 amol BSA digest.



MS/MS Fragment Analysis - AutoPSD

In the case of a questionable result from the digest, additional MS/MS (PSD) data can be collected from the same prepared sample. The new AutoPSD function facilitates the automated processing of this second analysis stage.




Go Proteineering!

The autoflex is part of the *PROTEINEER*TM suite – the most advanced hardware and software platform for protein identification and characterization.

Screening analysis - MALDI-TOF Peptide Mapping

The tryptic digest of a protein delivers peptides with masses of about 800 to 4,000 Da – the so-called “peptide map”, which is the basis for mass spectrometry-based protein identification. MALDI-TOF analysis of the digest and subsequent protein database search leads to a result list of proposed proteins.

The *PROTEINEER* suite supports this identification procedure as batch process. A central excel sheet contains all sample-specific parameters and controls all steps from sample preparation and MALDI-TOF analysis through to data interpretation.



The parent masses for PSD fragmentation are selected in the central Excel sheet.

biotools automatically annotates the sequence resulting from the search to the PSD spectrum. Thus the quality of the assignment is immediately visualized.

Production Sites

Bruker Daltonik GmbH

Bremen · Germany

Phone +49 (421) 22 05-200

Fax +49 (421) 22 05-103

E-Mail sales@bdal.de

Bruker Saxonia Analytik GmbH

Leipzig · Germany

Phone +49 (341) 24 31-30

Fax +49 (341) 24 31-404

E-Mail sales@bsax.de

Bruker Daltonics Inc.

Billerica · United States

Phone +1 (978) 663-36 60

Fax +1 (978) 667-59 93

E-Mail ms-sales@bdal.com

Sales and Service Centers

Europe

Belgium

Bruker Belgium S.A./N.V.

Bruxelles

Phone +32 (2) 726-76 26

Fax +32 (2) 726-82 82

E-Mail bruker@bruker.be

France

Bruker Daltonique S. A.

Wissembourg

Phone +33 (388) 73-68-00

Fax +33 (388) 73-68-79

E-Mail infomasse@bruker.fr

Great Britain

Bruker Daltonics Ltd.

Coventry

Phone +44 (24 76) 855-200

Fax +44 (24 76) 465-317

E-Mail sales@bruker.co.uk

Italy

Bruker Daltonics S.r.l

Milano

Phone +39 (02) 70 63-63 70

Fax +39 (02) 23 61-2 94

E-Mail bruker@bdal.it

Netherlands

Bruker Nederland NV

Wormer

Phone +31 (75) 628-52 51

Fax +31 (75) 628-97 71

E-Mail bruker@bruker.nl

Russia

Bruker Instruments

Moscow

Phone +7 (095) 137-67 51

Fax +7 (095) 137-67 51

E-Mail brukmos@cacr.ioc.ac.ru

Spain

Bruker Espanola S.A.

San Fernando de Henares

Phone +34 (91) 655 94 00

Fax +34 (91) 656 62 37

E-Mail bruker@bruker.es

Sweden

Bruker Daltonics Scandinavia AB

Täby

Phone +46 (8) 446-36 30

Fax +46 (8) 630-12 81

E-Mail ms@bruker.se

Switzerland

Bruker Daltonics AG

Fällanden

Phone +41 (1) 825-91 11

Fax +41 (1) 825-96 96

E-Mail daltonics@bruker.ch

America

Canada

Bruker Daltonics Ltd.

Milton, Ontario

Phone +1 (905) 876-46 41

Fax +1 (905) 876-44 21

E-Mail michael.mcdonell@bruker.ca

Mexico

Bruker Mexicana S.A. de C.V.

Mexico

Phone +52 (525) 630-57 47

Fax +52 (525) 630-57 46

United States

Bruker Daltonics Inc.

Fremont

Phone +1 (510) 683-43 00

Fax +1 (510) 490-65 86

E-Mail iso@bruker.com

Asia/Australia

Australia

Bruker Daltonics Pty. Ltd.

Alexandria

Phone +61 (2) 95 50-64 22

Fax +61 (2) 95 50-36 87

E-Mail pfb@bruker.com.au

China

Bruker Daltonics Inc.

Beijing

Phone +86 (010) 68 47-20 15

Fax +86 (010) 68 47-20 09

E-Mail hqwang@ihw.com.cn

India

Bruker India Scientific Pvt. Ltd.

Bombay

Phone +91 (22) 626-22 32

Fax +91 (22) 626-88 44

Japan

Nihon Bruker Daltonics K.K.

Tsukuba-Shi

Phone +81 (298) 52-12 34

Fax +81 (298) 58-03 22

E-Mail info@bruker.co.jp

Taiwan

Bruker Daltonics Inc.

Taipei

Phone +886 (2) 89 82-37 10

Fax +886 (2) 89 82-37 11

E-Mail ks.cheng@daltonics.bruker.com

Thailand

Bruker South East Asia

Bangkok

Phone +66 (2) 642-69 00

Fax +66 (2) 642-69 01

E-Mail bsea@bruker.com

Further representatives are available through the German or US headquarters.

www.bdal.de
www.bdal.com

BRUKER
DALTONICS®